

**Equipment Check List**  
**Adult and Pediatric Endotracheal Intubation Skill Station**

- \_\_\_ Adult training intubation manikin head
- \_\_\_ Pediatric training intubation manikin head
- \_\_\_ Latex gloves
- \_\_\_ Eye protection (goggles)
- \_\_\_ Assorted oropharyngeal airways of various sizes
- \_\_\_ Adult bag valve mask
- \_\_\_ Pediatric bag valve mask
- \_\_\_ Oxygen cylinder and tubing
- \_\_\_ Suction unit with soft and rigid suction catheters
- \_\_\_ Magill forceps (adult and pediatric)
- \_\_\_ Stethoscope
- \_\_\_ Laryngoscope handle and blades of various sizes and styles
- \_\_\_ Adult and pediatric endotracheal tubes of various sizes
- \_\_\_ Adult and pediatric stylets
- \_\_\_ Water soluble lubricant
- \_\_\_ 10cc syringe
- \_\_\_ One inch tape or endotracheal tube securing device
- \_\_\_ Stopwatch
- \_\_\_ Esophageal detector device (EDD)



	Points		Comments
10. Holds ET firmly in place until secured	1		
11. Correctly removes stylet and blade	1		
12. Inflates cuff with 5-10 cc of air prior to auscultating lung fields (uncuffed tubes are utilized for pediatric pts.)	1		
13. Confirms tube placement by using an esophageal detector device (EDD)	1		
14. Begins ventilation with BVM	1		
15. Checks tube placement by auscultating left and right lung fields, epigastrium & by observing chest rise	1*		
16. Secures endotracheal tube and inserts an oropharyngeal airway	1		
17. Re-verifies endotracheal tube placement as in steps 14 and 15	1		
18. If intubation attempt is >30 sec., the candidate ceases attempt & ventilates patient before reattempting intubation	1*		
<b>Total Points** -----&gt;</b>	<b>18</b>		

\* Failure to perform this step is a critical failure and requires the candidate to repeat the skill assessment station.

\*\* 15 points are required for a satisfactory rating. Failure to attain at least 15 points will require the candidate to repeat the skill assessment station.