

Equipment Check List
Automated External Defibrillator Skill Station

_____ Training manikin

_____ **Training** Automated External Defibrillator

_____ Adult bag-valve-mask

_____ Pocket face mask

Automated External Defibrillation Skill Assessment

Goals: The candidate will demonstrate proficiency using an automated external defibrillator and exhibit familiarity and comprehension of the Public Access Defibrillation protocols.

Candidate: _____ Preceptor: _____

Candidate's Recommended Renewal Date: _____

Status: SATISFACTORY
 REPEAT SKILL
 REFER TO REMEDIATION

	Points		Comments
1. Prepares equipment and exposes the patient's chest	1		
2. Attaches cables & positions pads appropriately	1		
3. Turns on defibrillator	1		
4. Ensures that no CPR is being done while machine analyzes.	1		
5. Confirms no carotid pulse	1		
6. Instructs all to "stand clear" and "clears" him/her self as indicated by the AED.	1*		
7. Rechecks for a pulse after the defibrillation.	1*		
8. Ensures that no CPR is being done while machine analyzes.	1		

	Points		Comments
9. Confirms no carotid pulse	1		
10. Instructs all to “stand clear” and “clears” him/her self as indicated by the AED.	1*		
11. Rechecks for a pulse after the defibrillation.	1		
12. Ensures that no CPR is being done while machine analyzes.	1		
13. Begins CPR when machine indicates to do so.	1		
Total Points** ----- >	13		

* Failure to perform this step is a critical failure and requires the candidate to repeat the skill assessment station.

** 13 points are required for a satisfactory rating. Failure to attain at least 13 points will require the candidate to repeat the skill assessment station.